Diagnostics is in our blood.

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

Profilaxis en HAG con/sin inhibidores

Emicizumab


FIGURE 2 Regulatory approval of emicizumab for inhibitor and non-inhibitor indications by year of approval. Data supplied courtesy of Roche Global, Switzerland

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



Table 1. Bleeding rates and proportion of patients with zero treated bleeds reported in clinical trials of emicizumab-kxwh.

Emicizumab-kxwh dose 1.5 mg/kg QW			
	HAVEN 1* [16]	HAVEN 2 [17]	HAVEN 3 [18]
Population	≥12 YO	<12 YO (and 12–17 YO if under 40 kg)	≥12 YO
Inhibitor status	Inhibitor	Inhibitor	Non-inhibitor
N patients	35	65	36
Treated bleeds			
ABR, model based (95% CI)	2.9 (1.69, 5.02)	0.3 (0.17, 0.50)	1.5 (0.9, 2.5)
% reduction (RR)#, P-value	87% (0.13), P < 0.0001	–	96%
Median ABR, calculated (IQR)	0.0 (0.0–3.7)	0.0 (0.0–0.0)	0.0 (0.0–2.5)
Zero bleeds, % (95% CI)	–	76.9% (64.8, 86.5)	56% (38, 72)
All bleeds			
ABR, model based (95% CI)	5.5 (3.58, 8.60)	3.2 (1.94, 5.22)	2.5 (1.6, 3.9)
% reduction (RR)#, P-value	80% (0.20), P < 0.0001	–	95%
Median ABR, calculated (IQR)	2.0 (0.0–9.9)	0.6 (0.00, 2.92)	0.6 (0.0–3.9)
Zero bleeds, % (95% CI)	–	49.2% (36.6, 61.9)	50% (33, 67)
Emicizumab-kxwh dose 3 mg/kg Q2W			
	HAVEN 2 [17]	HAVEN 3 [18]	HOHOEMI [19]
Population	<12 YO (and 12–17 YO if under 40 kg)	≥12 YO	<12 YO
Inhibitor status	Inhibitor	Non-inhibitor	Non-inhibitor
N patients	10	35	6
Treated bleeds			
ABR, model based (95% CI)	0.2 (0.03, 1.72)	1.3 (0.8, 2.3)	1.3 (0.59, 2.92)
% reduction (RR)#, P-value	–	97%	–
Median ABR, calculated (IQR)	0.0 (0.0–0.0)	0.0 (0.0–1.9)	1.4 (0.0–2.52)
Zero bleeds, % (95% CI)	90% (55.5, 99.7)	60% (42, 76)	33.3%
All bleeds			
ABR, model based (95% CI)	1.5 (0.62, 3.40)	2.6 (1.6, 4.3)	14.1 (7.63, 26.17)
% reduction (RR)#, P-value	–	94%	–
Median ABR, calculated (IQR)	0.0 (0.00, 2.81)	1.6 (0.0–4.0)	10.7 (5.51–20.37)
Zero bleeds, % (95% CI)	60% (26.2, 87.8)	40% (24, 58)	0
Emicizumab-kxwh dose 6 mg/kg Q4W			
	HAVEN 2 [17]	HAVEN 4 [14]	HOHOEMI [19]
Population	<12 YO (and 12–17 YO if under 40 kg)	≥12 YO	<12 YO
Inhibitor status	Inhibitor	Inhibitor/Non-inhibitor	Non-inhibitor
N patients	10	41 [§]	7
Treated bleeds			
ABR, model based (95% CI)	2.2 (0.69, 6.81)	2.4 (1.4, 4.3)	0.7 (0.18, 2.65)
Median ABR, calculated (IQR)	0.0 (0.0–3.26)	0.0 (0.0–2.1)	0.0 (0.0–1.73)
Zero bleeds, % (95% CI)	60% (26.2, 87.8)	56% (39.7, 71.5)	71.4%
All bleeds			
ABR, model based (95% CI)	3.8 (1.42, 10.11)	4.5 (3.1, 6.6)	21.8 (9.16, 51.80)
Median ABR, calculated (IQR)	1.6 (0.00, 4.84)	2.1 (0.0–5.9)	13.8 (4.58, 27.51)
Zero bleeds, % (95% CI)	50% (18.7, 81.3)	29% (16.1, 45.5)	14.3%

Mg, milligram; kg, kilogram; YO, year old; N, number; ABR, annualized bleeding rate; CI, confidence interval; RR, risk ratio; IQR, interquartile range.

* Bleeding rates of participants assigned to Group A (emicizumab-kxwh 1.5 mg/kg/wk) of HAVEN 1 who were receiving episodic treatment with BPAs prior to trial enrollment.

Percent risk reduction in bleeding compared to no emicizumab-kxwh prophylaxis.

§ 5 patients with inhibitors present at study entry.

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

Emicizumab

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
 <p style="text-align: center;">FVIIIa</p>	 <p style="text-align: center;">ACE910/Emicizumab</p>
Multiple sites of interaction	Single sites of interaction
High affinity for enzyme & substrate <i>(low to high nanomolar range)</i>	Low affinity for enzyme & substrate <i>(micromolar range)</i>
Specific for FIXa and FX <i>(no binding to FIX and FXa)</i>	No distinction between zymogen and enzyme <i>(FIX vs FIXa and FX vs FXa)</i>
Full cofactor activity <ul style="list-style-type: none"> - <i>promotes phospholipid binding</i> - <i>stabilizes FIXa active site</i> - <i>bridges FIXa to FX</i> 	Partial cofactor activity <ul style="list-style-type: none"> - <i>bridges FIXa to FX</i>
Enzyme and substrate are in excess over cofactor	Antibody is in excess over enzyme and substrate
FVIIIa has on/off mechanism	Emicizumab has no on/off mechanism
High level of self-regulation	Low level of self-regulation

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

Test results affected by emicizumab	Test results unaffected by emicizumab
APTT reduced	
PT reduced [described as a clinically insignificant effect]	PT-based assays
APTT-based single factor assays (increased)	Thrombin time
APTT-based inhibitor screens (false negative)	Clauss fibrinogen
APTT-based factor inhibitor titres (false negative)	Immuno-based assays
Chromogenic FVIII assays using human-derived FIXa and FX (increased)	Chromogenic FVIII assays using bovine-derived FIXa and FX
Activated clotting time (ACT) (possibly reduced)	Chromogenic assays for coagulation factors or treatments other than FVIII.

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
FVIII activity assays	<u>FVIII OSA (aPTT-based)</u>	FVIII: sensitive Emicizumab: oversensitive; normalizes at very low concentrations of emicizumab ^{16,36}	Cannot be used to measure emicizumab or FVIII activity
	<u>Modified FVIII OSA calibrated against emicizumab</u>	FVIII: slightly sensitive Emicizumab: sensitive ^{41,46}	<u>May be used to measure emicizumab plasma levels</u>
	<u>FVIII CSA (bovine components: Siemens, etc.)</u>	FVIII: sensitive Emicizumab: insensitive; emicizumab does not accelerate FXa formation by bovine FIXa ⁵⁷	Cannot be used to measure emicizumab activity. <u>May be used to measure FVIII activity without interference from emicizumab in the sample</u>
	<u>FVIII CSA (human components: Hyphen BIOPHEN)</u>	FVIII: sensitive Emicizumab: sensitive (within the dynamic range of the assay) ⁴¹	Reported <u>FVIII cannot be viewed as equivalent to FVIII activity measured in patients treated with FVIII. May only provide a relative indication of the procoagulant activity of emicizumab</u>

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Inhibitor assays	Bethesda assays using FVIII OSA	FVIII: sensitive Emicizumab: false negative; as emicizumab is not deactivated by heat, it drives coagulation via human FIX and FX in human plasma, regardless of presence of inhibitors to FVIII ³¹	<u>Cannot be used to measure FVIII inhibitor titers in the presence of emicizumab</u>
	CBA	FVIII: sensitive Emicizumab: insensitive; this assay uses a bovine protein-based chromogenic FVIII assay, which is not affected by emicizumab ⁴²	<u>Can be used to measure FVIII inhibitor titers in the presence of emicizumab</u>



Abbreviations: aPTT, activated partial thromboplastin time; CBA, chromogenic Bethesda Assay; CSA, chromogenic assay; F, coagulation factor; OSA, one-stage clotting assay.

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



	Standard assay	Sensitivity to FVIII/emicizumab	Recommendations for use in the presence of emicizumab
Global coagulation assays	<u>aPTT</u>	FVIII: sensitive Emicizumab: oversensitive; normalizes at very low concentrations of emicizumab ³⁹	Cannot be used to measure the hemostatic potential of a patient
	<u>ROTEM</u>	FVIII: sensitive Emicizumab: sensitive; EXTEM has little-to-no response in the presence of emicizumab and INTEM shows marked response in the presence of emicizumab, but no dose-dependent effect, while NATEM shows a significant dose-dependent response to emicizumab ⁵²	May be used to measure global coagulation in the presence of emicizumab
	<u>TGA</u>	FVIII: sensitive Emicizumab: sensitive; shows a linear, dose-dependent increase in thrombin generation when triggered with FXIa ⁴⁸	May be used to measure global coagulation in the presence of emicizumab

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



N = 58 muestras

Estudio HAG con/sin inhibidor



Niveles de emicizumab vs Lag Time y Time to Peak

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